REvision

AZURE HT Technology

EN

endo★star

Endostar REvision

Instruction for use

Endostar REvision

Endostar REvision made with AZURE HT Technology by Poldent - Innovative heat treatment technology created by Poldent.

1. Important information regarding the system

The Endostar REvisions is a rotary instrument for efficient removal of old root canal fillings during a revision endodontic procedure with the crown down technique. The set contains three files marked with one (#1), two (#2) or three (#3) stripes on the shaft (sizes: 08/30, 06/25, 04/20 respectively).

Endostar REvision files are manufactured from the highest quality nickel-titanium alloy, which additionally underwent a special heat treatment called AZURE HT Technology by Poldent, thus achieving very high flexibility and durability.

> The handpiece, which can provide 300 rpm, should be used. The operating speed of the handpiece should be constant throughout the shaping process.

> Do not apply excessive force. An up-and-down motion should be used when operating the files.

> Shaping time should be as short as possible.

> Always use a lubricating agent when shaping the canal.

> The files are very sharp and should be used very carefully, with little force and without excessive "pushing" down the canal.

> Operate the instruments and handpieces according to their operating instructions (especially torque and speed settings).

> Use the type and amount of instruments that is truly needed in a given clinical situation.

> Before using the instruments, be sure to see them working outside the oral cavity to check for deformations and/or cracks.

> Dispose of as medical waste.

> The Endostar REvision files may remain bend and may not straighten at room temperature as the non-modified NiTi alloy files do. This is a normal feature of the instrument.

> The Endostar REvision files can be pre-bent, the same way as the steel files, before inserting them into the root canal in order to bypass the existing ledges.

> It is also acceptable to insert pre-bent file into the canal and then start the micromotor, this simplifies the access to the root canals in molars.

2. Recommended movements

All instruments has been designed and manufactured in such a way that it can be used in three types of movements depending on the individual preferences of the dentist, the case diagnosis and the type of a handpiece available in the dental practice.

> **Rotary movement** - the instrument rotates continuously 360° in a clockwise direction (CW - Clock Wise).

> **Reciprocal right cutting movement** - the instrument performs alternating movements: clockwise (CW) and anti-clockwise (CCW - Counter Clock Wise) except that the CW movement angle must be larger than CCW angle, e.g. 90° CW and 30° CCW.

It is recommended that the rotation in the CW direction should be between of 90° to 270° and in the CCW direction between 30° to 90°, so that the net rotation in the CW direction in each cycle is between 60° to 240°, that means a full 360° CW rotation is achieved after 1.5 to 6 cycles.

> **Complex movement** - it is a kind of movement that combines the rotary movement with the reciprocal movement. After inserting file into the root canal, the file performs a rotary motion, and if the resistance for the file in the canal is too high, the rotary motion changes to the reciprocating movement. When the resistance decreases, the rotary motion returns. An example of this is an OTR movement.

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